



CUSTOMED



It's all about healing

Sterilization Wrappers

Our sterilization wrappers come in a variety of sizes and grades to meet your specific needs. They are available with the following specifications: 100% Polypropylene, 100% SMS, 52 gsm, anti-static and anti-repellent.



Size	Item No.:	Cat. No.	Packaging:
10" x 10"	712057470010	900078	2000/CS
15" x 15"	712057470015	900079	1000/CS
20" x 20"	712057470020	900081	1000/CS
24" x 24"	712057470024	900082	750/CS

SMS Latex Free 100% Polypropylene

Size	Item No.:	Cat. No.	Packaging:
30" x 30"	712057470030	900083	250/CS

SMS Latex Free 100% Polypropylene

Size	Item No.:	Cat. No.	Packaging:
40" x 40"	712057470040	900084	250/CS
45" x 45"	712057470044	900085	250/CS

SMS Latex Free 100% Polypropylene

CUSTOMED



It's all about healing

10805 Southport Drive, Unit #1
Orlando, FL 32824
1.800.981.2679 | 1.407.850.5558
www.customedhealing.com